

# 637. Thiamethoxam

## HPLC method

### **CIPAC Collaborative Trial according to CIPAC Information Sheet N° 293**


Dr. Sven Adolph  
Syngenta Crop Protection Münchwilen AG  
CH-4333 Münchwilen

Switzerland

May 2012

## THIAMETHOXAM

### 637/TC/M/-



<i>ISO Common Name:</i>	Thiamethoxam
<i>Chemical Name:</i>	3-(2-chloro-thiazol-5-ylmethyl)-5-methyl-[1,3,5]oxadiazinan-4-ylidene-N-nitro-amine
<i>CAS-Number:</i>	153719-23-4
<i>Molecular mass:</i>	291.7
<i>Empirical formula:</i>	C <sub>8</sub> H <sub>10</sub> ClN <sub>5</sub> O <sub>3</sub> S
<i>m.p.</i>	139.1°C
<i>b.p.</i>	Decomposition occurs at about 147 °C before boiling point is reached

**1 Sampling.** Take at least 100 g.

### 2 Identity test

**2.1 HPLC.** Use the reversed phase HPLC method described below in section 3.1. The relative retention time of the Thiamethoxam peak in the sample solution should not deviate by more than 2% from that of the calibration solution.

**2.2 GC.** Use the capillary GC method below. The relative retention time of the Thiamethoxam peak in the sample solution should not deviate by more than 2% from that of the calibration solution.

#### REAGENT

5ml water / 95ml acetone

#### APPARATUS

*Gas chromatography* equipped with flame ionisation detector.

Capillary column fused silica, length 30 m x 0.25 (i.d.) mm, film thickness: 0.5 µm, coated with crosslinked (5%-Phenyl)-methylpolysiloxane (DB-5 or equivalent).

Data system

## PROCEDURE

(a) Gas chromatographic conditions (typical):

<i>Column</i>	Fused silica, length 30 m x 0.25 mm internal diameter. DB-5 or equivalent phase, with film thickness of 0.5 µm crosslinked (5%-Phenyl)-methylpolysiloxane
<i>Injection system</i>	
Injector	split injection
Injection volume	1 µl
Split ratio	1:20
<i>Detector</i>	flame ionisation
<i>Temperatures</i>	
Injection port	260 °C
Detector	300 °C
Oven programme	temp 1 60°C, hold 1 min, ramp rate 20°C/min temp 2 170°C, hold 0 min, ramp rate 3°C/min temp 3 200°C, hold 0 min, ramp rate 20°C/min temp 4 300°C, hold 5 min,
<i>Gas flow rates</i>	
<i>Column:</i>	
Hydrogen (carrier)	2 ml/min (typically 11.7 psi at 60°C), run at constant flow
<i>Detector:</i>	
Air	300 ml/min
Hydrogen	30 ml/min
Nitrogen (make up)	30 ml/min
<i>Retention times</i>	
Thiamethoxam	18.3 min (approximate)

*Sample preparation:* Weigh (to the nearest 0.1 mg) sufficient sample (*w* mg) to contain about 200 mg Thiamethoxam into a volumetric flask (100 ml). Add acetone (about 50 ml) and water (5 ml) and place the flask in an ultrasonic bath for 5 min. Allow to cool to ambient temperature and fill to the mark with acetone. Mix thoroughly.

*Calibration solution:* Weigh (to the nearest 0.1 mg) sufficient sample (*w* mg) to contain about 200 mg Thiamethoxam into a volumetric flask (100 ml). Add acetone (about 50 ml) and water (5 ml) and place the flask in an ultrasonic bath for 5 min. Allow to cool to ambient temperature and fill to the mark with acetone. Mix thoroughly

### 3 Thiamethoxam

#### OUTLINE OF METHOD

Thiamethoxam content is determined (g/kg) using reversed phase high performance liquid chromatography incorporating UV detection at 230 nm with an external standard calibration.

#### 3.1 Determination of Thiamethoxam by reversed phase HPLC

##### REAGENTS

*Thiamethoxam* reference standard of known content

*Water* HPLC grade

*Acetonitrile* HPLC grade

*Phosphoric acid (85%)* analytical grade

*Eluent A:* 0.1% v/v aqueous phosphoric acid

*Eluent B:* acetonitrile

*Diluent:* 0.1% v/v aqueous phosphoric acid / acetonitrile 1 + 4 (v/v)

*Calibration solution.* Weigh in duplicate (to the nearest 0.1 mg) 45 - 55 mg of the Thiamethoxam reference standard (s mg) into separate volumetric flasks (50 ml). Add about 35 ml acetonitrile and 10 ml 0.1% v/v aqueous phosphoric acid and place the flasks in an ultrasonic bath for 5 minutes. Allow the solutions to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Transfer 25.0 ml of this solution in a 100 ml volumetric flask, add diluent (0.1% v/v aqueous phosphoric acid / acetonitrile, in a ratio of 1/4 (v/v)) and make up to the volume at room temperature (calibration solutions C<sub>1</sub> and C<sub>2</sub>).

##### APPARATUS

*High performance liquid chromatograph* equipped with a detector suitable for operation at 230 nm (UV-detection) and an injection system capable to inject 5 µl.

*Liquid chromatographic column* stainless steel, 75 x 4.6 mm (i.d), packed with Nucleodur C<sub>18</sub> ec (3 µm) Macherey-Nagel or equivalent with the same selectivity.

*Electronic integrator or data system*

*Ultrasonic bath*

## PROCEDURE

(a) *Chromatographic conditions* (typical):

<i>Column temperature</i>	40°C
<i>Flow rate</i>	1.5 ml/min
<i>Detector wavelength</i>	230 nm
<i>Injection volume</i>	5 µl
<i>Retention time</i>	approximately 2.6 minutes

<i>Gradient program:</i>	time [minutes]	0.1 % aqueous phosphoric acid [%]	acetonitrile [%]
	0	90	10
	4	30	70
	10	5	95
	11	5	95
	11.1	90	10
	14	90	10

(b) *Equilibration of the system.* Pump sufficient mobile phase through the column to equilibrate the system. Inject 5 µl portions of the calibration solution C<sub>1</sub> and repeat the injections until retention times and peak areas deviate by less than ± 1 % from the mean for three successive injections.

(c) *Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) 45 - 55 mg of Thiamethoxam tech. (w mg) into a volumetric flask (50 ml). Add about 35 ml acetonitrile and 10 ml 0.1 % v/v aqueous phosphoric acid and place the flasks in an ultrasonic bath for 5 minutes. Allow the solutions to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Transfer 25.0 ml of this solution in a 100 ml volumetric flask, add diluent (0.1% v/v aqueous phosphoric acid / acetonitrile, in a ratio of 1/4 (v/v)) and make up to the volume at room temperature (sample solutions S<sub>1</sub> and S<sub>2</sub>).

(d) *Determination.* Inject 5 µl portions of the second calibration solution (C<sub>2</sub>) for two successive injections. The mean response factor for this solution should deviate by no more than 1% from those for the first calibration solution (C<sub>1</sub>) (see paragraph (b) *Equilibration of the system*), otherwise the calibration solutions should be prepared again.

Inject in duplicate 5 µl portions of each sample solution (S<sub>1</sub>, S<sub>2</sub>, ..., etc.) bracketing them by duplicate injections of the calibration solution (C<sub>1</sub>) using the following sequence:

C<sub>1</sub>, C<sub>1</sub>, S<sub>1</sub>, S<sub>1</sub>, S<sub>2</sub>, S<sub>2</sub>, C<sub>1</sub>, C<sub>1</sub>...

Determine the peak area of Thiamethoxam.

(e) *Calculation.* Calculate the response factors (*f*) from the calibration solutions bracketing the injections of the sample solutions. Average the response factors of the calibration solutions preceding and following the sample solution injections. These must agree within ± 1 % of the average otherwise repeat the determination. Calculate the content of the sample solutions.

$$f_i = \frac{s \times P}{H_s}$$

$$\text{Thiamethoxam content} = \frac{H_w \times f}{w} \text{ [g/kg] } (M)$$

where:

*f<sub>i</sub>* = single response factor

*f* = average response factor

*H<sub>s</sub>* = peak area of Thiamethoxam in the calibration solution

*H<sub>w</sub>* = peak area of Thiamethoxam in the sample solution

*s* = mass of Thiamethoxam reference standard in the calibration solution (mg)

*w* = mass of sample taken (mg)

*P* = purity of Thiamethoxam reference standard (g/kg)

**THIAMETHOXAM WATER DISPERSIBLE GRANULE**  
637/WG/M/-

**1 Sampling.** Take at least 100 g.

**2 Identity test**

**2.1 HPLC.** As for Thiamethoxam technical 637/TC/M/2.1

**2.2 GC** As for Thiamethoxam technical 637/TC/M/2.2

**3 Thiamethoxam**

Same approach as for Thiamethoxam technical 637/TC/M/3

**3.1 Determination of Thiamethoxam by reversed phase HPLC**

As for Thiamethoxam technical 637/TC/M/3.1 except

*Sample filtering device* with a membrane filtration unit compatible with organic solvents and a 0.45 µm pore diameter.

(c) *Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample to contain 45 - 55 mg Thiamethoxam (w mg, equal to 180 – 220 mg Thiamethoxam formulation WG 25) into a volumetric flask (50 ml). Suspend with 10 ml 0.1% v/v aqueous phosphoric acid. Add about 35 ml acetonitrile and place the flasks in an ultrasonic bath for 5 minutes. Allow the solutions to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Filter solutions through a 0.45 µm filter. Transfer 25.0 ml of this solution in a 100 ml volumetric flask, add diluent (0.1% v/v aqueous phosphoric acid / acetonitrile, in a ration of 1/4 (v/v)) and make up to the volume at room temperature (sample solutions S<sub>1</sub> and S<sub>2</sub>).

## 4 Suspensibility

REAGENTS AND APPARATUS as for 637/TC/M/- and MT 184.

add at:

APPARATUS

*Sample filtering device* with a membrane filtration unit compatible with organic solvents and a 0.45 µm pore diameter

PROCEDURE

(a) *Preparation of suspension and determination of sedimentation.* MT 184.

(b) *Determination of Thiamethoxam in the bottom 25 ml of suspension.* After removal of the top 225 ml of suspension, transfer the remaining 25 ml quantitatively into a volumetric flask (100 ml) with acetonitrile (about 50 ml). Place the flask in an ultrasonic bath for 5 minutes. Allow to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Dilute 10.0 ml of the solution with the eluent to a volume of 100 ml. Mix thoroughly. Clear the suspension by filtration through a 0.45 µm filter prior to injection. Determine the mass (Q g) of Thiamethoxam according to 637/WG/M/3.1, using a calibration solution with the appropriate final concentration of Thiamethoxam standard.

(c) *Calculation*

$$\text{Suspensibility} = \frac{111(c - Q)}{c} \%$$

where:

c = mass of Thiamethoxam in the sample taken for the preparation of the suspension (g)

Q = mass of Thiamethoxam in the bottom 25 ml of suspension (g)



**THIAMETHOXAM FLOWABLE CONCENTRATE**  
637/SC/M/-

**1 Sampling.** Take at least 100 g.

**2 Identity test**

**2.1 HPLC.** As for Thiamethoxam technical 637/TC/M/2.1

**2.2 GC** As for Thiamethoxam technical 637/TC/M/2.2

**3 Thiamethoxam**

Same approach as for Thiamethoxam technical 637/TC/M/3

**3.1 Determination of Thiamethoxam by reversed phase HPLC**

As for Thiamethoxam technical 637/TC/M/3.1 except

*Sample filtering device* with a membrane filtration unit compatible with organic solvents and a 0.45 µm pore diameter.

(c) *Sample preparation.* Homogenize sample thoroughly. Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample to contain 45 - 55 mg Thiamethoxam (w mg, equal to 210 – 250 mg Thiamethoxam formulation SC 240) into a volumetric flask (50 ml). Suspend with 10 ml 0.1% v/v aqueous phosphoric acid. Add about 35 ml acetonitrile and place the flasks in an ultrasonic bath for 5 minutes. Allow the solutions to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Filter solutions through a 0.45 µm filter. Transfer 25.0 ml of this solution in a 100 ml volumetric flask, add diluent (0.1% v/v aqueous phosphoric acid / acetonitrile in a ratio of 1/4 (v/v)) and make up to the volume at room temperature (sample solutions S<sub>1</sub> and S<sub>2</sub>).

---

**THIAMETHOXAM FLOWABLE CONCENTRATE FOR SEED TREATMENT**  
637/FS/M/-

**1 Sampling.** Take at least 100 g.

**2 Identity test**

**2.1 HPLC.** As for Thiamethoxam technical 637/TC/M/2.1

**2.2 GC** As for Thiamethoxam technical 637/TC/M/2.2

**3 Thiamethoxam**


Same approach as for Thiamethoxam technical 637/TC/M/3

**3.1 Determination of Thiamethoxam by reversed phase HPLC**


As for Thiamethoxam technical 637/TC/M/3.1 except

*Sample filtering device* with a membrane filtration unit compatible with organic solvents and a 0.45 µm pore diameter.


*(c) Sample preparation.* Homogenize sample thoroughly. Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample to contain 45 - 55 mg Thiamethoxam (w mg, equal to 130 – 160 mg Thiamethoxam formulation FS 350) into a volumetric flask (50 ml). Suspend with 10 ml 0.1% v/v aqueous phosphoric acid. Add about 35 ml acetonitrile and place the flasks in an ultrasonic bath for 5 minutes. Allow the solutions to cool to ambient temperature and fill to the mark with acetonitrile. Mix thoroughly. Filter solutions through a 0.45 µm filter. Transfer 25.0 ml of this solution in a 100 ml volumetric flask, add diluent (0.1% v/v aqueous phosphoric acid / acetonitrile in a ratio of 1/4 (v/v)) and make up to the volume at room temperature (sample solutions S<sub>1</sub> and S<sub>2</sub>).




**Fig 1** Typical HPLC-chromatogram of Thiamethoxam technical material




**Fig 2** Typical HPLC-chromatogram of Thiamethoxam formulation WG 25



**Fig 3** Typical HPLC-chromatogram of Thiamethoxam in formulation SC 240



**Fig 4** Typical HPLC-chromatogram of Thiamethoxam in formulation FS 350



**Fig 5** Typical GC-chromatogram of Thiamethoxam technical material