



Jiangsu Yangnong Chemical Co., Ltd.
NO.39 WenfengRoad.
Yangzhou
Jiangsu
CHINA

Tel: +86-514-85889958
Fax: +86-514-85889900
e-mail : trade@yangnong.net

5040m. d-Tetramethrin

GC and LC method 5040m


CIPAC Small Scale Collaborative Trial

Jiangsu Yangnong Chemical Co., Ltd.
YANGZHOU

China

June 2016

D-TETRAMETHRIN
5040/TC/M



ISO common name: Not available

Other names: d-tetramethrin

Chemical name: (1,3,4,5,6,7-Hexahydro-1,3-dioxo-2H-isoindol-2-yl)methyl (1R-trans)-2,2-dimethyl-3-(2-methylprop-1-enyl)cyclopropanecarboxylate

CAS-Number: 1166-46-7

RMM: 331.42

Empirical formula: C₁₉H₂₅NO₄

m.p. 48 °C

V.p. 9.5x 10⁻⁴ Pa at 30°C

Solubility
In water at 20 °C: insoluble
In toluene: >200 g/l at 20°C
In acetone: >200 g/l at 20°C

Stability Stable at ambient temperature

Description Yellow to brown liquid

Note: d-tetramethrin is a mixture of the isomers (1R-trans, R), (1R-trans, S), (1R-cis, R) and (1R-cis, S) of tetramethrin in an approximate ratio of 4:4:1:1. In practice the trans isomer range is 75-85 % and the cis isomer range is 15-25 %.

D-TETRAMETHRIN TECHNICAL
5040m/TC/M/-

1 Sampling. Take at least 100 g.

2 Identity test

2.1 Infrared. Prepare a film between KBr plates and scan from 4000 to 400 cm⁻¹. The spectrum produced from the sample should not differ significantly from that of the standard.

2.2 HPLC. Use the normal phase HPLC method below. The relative retention time of the respective peak in the sample solution should not deviate by more than 1.5% from that of the calibration solution.

n-hexane HPLC grade

Ethanol HPLC grade

Diethylamine AR grade

Mobile phase. Add 1 ml diethylamine and 70 ml ethanol to 930 ml *n*-hexane and mix. Mixing of the mobile phase may also be accomplished using a binary solvent HPLC pump. Filter and degas the mobile phase prior to use.

APPARATUS

High performance liquid chromatograph equipped with a normal-flow pump, normal-temperature chiral column compartment, a sample injector capable of injecting 3 µl aliquots, a UV detector (230 nm) and digital integrator or other data-handling capability. A chilled autosampler should be used, if available (but is not required).

HPLC column CHIRALCEL[®] AY-H, 250 × 4.6 mm (i.d.), 5 µm particle size;

Ultrasonic bath

(a) *Operating conditions* (typical)

Column 250 × 4.6 mm (i.d.), packed with CHIRALCEL[®] AY-H, 5 µm particle size

Mobile phase *n*-hexane-ethanol-diethylamine, 930+70+1 (v/v)

Column temperature 35 °C

Flow rate 1.0 ml/min

Injection volume 3 µl

Detector wavelength 230 nm

Run time approximately 20 min

Retention time 1*S*-cis isomer: about 7.7 min

1*R*-cis isomer: about 9.0 min

1*R*-trans isomer: about 10.0 min

1*S*-trans isomer: about 11.4 min

(b) *Sample preparation.* Melt the sample and homogenize thoroughly by stirring prior to sampling. Prepare a solution of approximately 250 mg/l of the sample in mobile phase. Inject 3 µl of this solution.

(c) *Equilibration of the system.* Inject 3 µl portions of the sample solution and repeat the injections until retention times and peak areas vary by less than 1.0% of the mean for 3 successive injections.

(d) *Determination.* Inject 3 µl of the sample solution and determine the respective peak areas.

(e) *Calculation.* Calculate the trans isomer and 1*R* isomer fraction percentage of the four isomers.

$$\text{Trans isomer fraction percentage} = \frac{H_c + H_d}{H_a + H_b + H_c + H_d} \times 100\%$$

$$\text{1R isomer fraction percentage} = \frac{H_b + H_c}{H_a + H_b + H_c + H_d} \times 100\%$$

Where:

H_a = area of the 1S-cis isomer

H_b = area of the 1R-cis isomer

H_c = area of the 1R-trans isomer

H_d = area of the 1S-trans isomer

3 d-Tetramethrin

OUTLINE OF METHOD. d-Tetramethrin is determined by gas chromatography with internal standardization and flame ionization detection.

REAGENTS

Dichloromethane

Tetramethrin working standard with known content

Di-(2-ethylhexyl)-O-phthalate (DOP) internal standard. Must not show peaks with the same retention times as Tetramethrin.

Internal standard solution. Dissolve DOP (1.0 g) in dichloromethane (100 ml). Ensure that a sufficient quantity of this solution is prepared for all samples and calibration standards to be analysed.

Calibration solution. Prepare calibration solutions in duplicate. Weigh (to the nearest 0.1 mg) 90 to 110 mg (s mg) of tetramethrin working standard into a volumetric flask (20 ml). Add by pipette internal standard solution (5.0 ml) and dissolve completely. Make up to volume with dichloromethane and mix well (solutions C_A and C_B).

APPARATUS

Gas chromatograph capable of operating in the range 240 °C equipped with a flame ionisation detector, a split injector and an autosampler

Column quartz, 30 m × 0.25 mm (i.d.), coated with 50% (trifluoropropyl)-methylpolysiloxane (e.g DB-210), film thickness 0.25 µm

Electronic integrator or data system

PROCEDURE

(a) *Chromatographic conditions* (typical)

Column quartz, 30 m × 0.25 mm (i.d), coated with 50% (trifluoropropyl) - methylpolysiloxane (e.g DB-210), film thickness 0.25 µm

Injection System

Injector Split injection

Split ratio 30:1

Detector Flame ionization

Temperatures

Injector 250 °C

Detector 280 °C

Oven 240 °C

Injection volume 1 µl

Gas flow rates

d-Tetramethrin 5040m

Helium (carrier)	1.6 ml/min
Hydrogen	about 40 ml/min
Air	about 400 ml/min
Helium (make up)	about 30 ml/min
<i>Run time</i>	about 10 min
<i>Retention times</i>	d- tetramethrin: about 7.4 min + 8.0 min
DOP	about 5.7 min

- (b) *Linearity check.* Check the linearity of the detector response by injecting 1 μ l of solutions with tetramethrin concentrations 0.5, 1 and 2 times that of the calibration solution before conducting the analysis.
- (c) *Equilibration of the system.* Prepare two calibration solutions. Inject 1 μ l portions of the first calibration solution C_A and repeat the injections until the peak area ratio of tetramethrin to the internal standard obtained from two consecutive injections of C_A deviate by less than 1.0%. Then the ratio for the second calibration solution C_B should not deviate by more than 1.0% from that for C_A . Otherwise prepare new calibration solutions.
- (d) *Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) 90 to 110 mg d- tetramethrin (w mg) of sample into a volumetric flask (20 ml). Add by pipette internal standard solution (5.0 ml) and dissolve completely. Make up to volume with dichloromethane and mix well (solutions S_A and S_B).
- (e) *Determination.* Injection in duplicate 1 μ l portions of each sample solution bracketing them by injections of the calibration solutions as follows: calibration solution C_A , sample solution S_A , sample solution S_A , calibration solution C_B , sample solution S_B , sample solution S_B , calibration solution C_A , and so on. Measure the relevant peak areas.
- (f) *Calculation.* Calculate the mean value of each pair of response factors bracketing the two injections of a sample and use this value for calculating the d- tetramethrin contents of the bracketed sample injections.

$$f_i = \frac{I_r \times s \times P}{H_s}$$

$$\text{Content of d- tetramethrin} = \frac{f \times H_w}{I_q \times w} \text{ g / kg}$$

Where:

f_i = individual response factor

f = mean response factor

H_s = peak area of d- tetramethrin in the calibration solution

H_w = peak area of d- tetramethrin in the sample solution

I_r = peak area of the DOP in the calibration solution

I_q = peak area of the DOP in the sample solution

s = mass of tetramethrin in the calibration solution (mg)

w = mass of sample taken (mg)

P = purity of tetramethrin standard (g/kg)

d-Tetramethrin 5040m




Figure 1 Gas Chromatogram of d- tetramethrin technical




Figure 2 Infrared spectrum of d- tetramethrin




Figure 3 HPLC Chromatogram of tetramethrin technical




Figure 4 HPLC Chromatogram of d-tetramethrin technical