

# **Flupyradifurone**

## **HPLC Method**

**CIPAC Collaboration Trial  
according to  
CIPAC Information Sheet No 308**

by


Alexandra Michel  
Crop Science Division  
Bayer Aktiengesellschaft  
Alfred-Nobel-Str. 50, Building 6820  
40789 Monheim am Rhein  
Germany

May 2017

**FLUPYRADIFURONE**  
**987**

Chemical name: 4-[[[(6-chloropyridin-3-yl)methyl](2,2-difluoroethyl)amino]furan-2(5H)-on  
 ISO common name: Flupyradifurone  
 CAS-No.: 951659-40-8

Structure:



Molecular mass: 288.7  
 Empirical formula: C<sub>12</sub> H<sub>11</sub> Cl F<sub>2</sub> N<sub>2</sub> O<sub>2</sub>

m.p.: 69 °C  
 b.p.: Decomposition occurs at about 270°C before boiling point is reached  
 Solubility: 1,2-dichloromethane, ethylacetate, methanol > 250 g/L; toluene 3.7 g/L (all at 20°C)

Description: The technical product is a white to beige powder.

**FLUPYRADIFURONE TECHNICAL**  
**987/TC/M/-**

**1 Sampling.** Take at least 100 g\*. Grind the sample thoroughly in a mortar.  
 \*(for this trial less amount is provided; please grind the entire sample)

**2 Identity tests**

**2.1 HPLC.** Use the HPLC method described below. The relative retention time of Flupyradifurone in the sample solution should not deviate by more than 2% from that of the calibration solution.

**2.2 UV spectrometry.** Record the UV spectrum during the HPLC determination. The UV spectrum obtained from the sample should not differ significantly from that of the standard. (Fig. 1)

**2.3 Infrared.** Prepare by direct application pure flupyradifurone and the sample on an ATR unit. Scan from 4000 to 550 cm<sup>-1</sup>. The spectrum produced from the sample should not differ significantly from that of the standard. (Fig. 2)

### 3 Flupyradifurone

#### OUTLINE OF THE METHOD.

Flupyradifurone content is determined (g/kg) by reversed phase high performance liquid chromatography using UV detection at 280 nm and external standard calibration.

#### 3.1 Determination of Flupyradifurone by reversed phase HPLC

##### REAGENTS

Flupyradifurone reference standard of known content  
 Acetonitrile (HPLC grade)  
 Phosphoric acid 85 % (puriss. p. a.)  
 Purified water (HPLC grade)  
 Eluent A: 10 mMol phosphoric acid in 1 L purified water  
 Eluent B: acetonitrile

##### APPARATUS

High performance liquid chromatograph equipped with an injection system capable to inject 5  $\mu$ L and an UV spectrophotometric detector operated at 280 nm.  
 Liquid chromatography column, stainless steel, 50 x 4.6 (i.d.) mm, packed with Kinetex C 18; 2.6  $\mu$ m or equivalent with the same selectivity.  
 Electronic integrator or data system  
 Ultrasonic bath

##### PROCEDURE

###### (a) Operating conditions (typical):

Flow rate: 2 mL/min  
 Column temperature: 40°C  
 Injection volume: 5  $\mu$ L  
 Detector wavelength: 280 nm

Mobile phase:

Time [min]	10 mMol phosphoric acid in 1 L purified water [%v/v]	Acetonitrile [%v/v]
0.0	90	10
5.0	85	15
5.1	05	95
6.5	05	95
6.6	90	10
8.5	90	10

Retention time: approximately 4.1 min

(b) *Equilibration of the system.* Pump sufficient mobile phase through the column to equilibrate the system. Inject 5  $\mu$ L portions of the calibration solution C1 and repeat the injections until retention times and peak areas deviate by less than  $\pm 1$  % from the mean for three successive injections.

(c) *Calibration solution.* Weigh in duplicate (to the nearest 0.1 mg) approximately 50 mg (*s* mg) of the Flupyradifurone reference standard into separate volumetric flasks (100 mL). Add 50 mL acetonitrile and place the flasks in an ultrasonic bath for 15 min. Make up the flasks with purified water to just below the calibration mark and allow to cool to ambient temperature. Fill to the mark with purified water and mix thoroughly (Calibration solutions C1, C2) (Fig. 1).

(d) *Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample (*w* mg) to contain about approximately 50 mg (*w* mg) of Flupyradifurone into separate volumetric flasks (100 mL). Add 50 mL acetonitrile and place the flasks in an ultrasonic bath for 15 min. Make up the flasks with purified water to just below the calibration mark and allow to cool to ambient temperature. Fill to the mark with purified water and mix thoroughly (Sample solutions S1, S2) (Fig. 2).

(e) *Determination.* Inject in duplicate each sample solution and bracket a series of sample solution injections by injections of the calibration solutions as follows: calibration solution 1, calibration solution 2, calibration solution 1, sample solution 1, sample solution 1, sample solution 2, sample solution 2, calibration solution 1, ... (C1, C2, C1, S1, S1, S2, S2, C1, ...).

Determine the peak areas of Flupyradifurone.

(f) *Calculation*

Calculate the response factors from the calibration solutions bracketing the injections of the sample solutions. Average the response factors of the calibration solutions preceding and following the sample solution injections. These must agree within  $\pm 1\%$  of the average otherwise repeat the determination. Calculate the content of the sample solutions.

$$f_i = \frac{s \times P}{H_s}$$

$$\text{Flupyradifurone content (g/kg)} = \frac{H_w \times f}{w}$$

Where:

$f_i$  = single response factor

$f$  = average response factor

$H_s$  = peak area of Flupyradifurone standard in the calibration solution

$H_w$  = peak area of Flupyradifurone in the sample solution

$s$  = weight of the Flupyradifurone standard in the calibration solution (mg)

$w$  = weight of the sample (mg)

$P$  = purity of the Flupyradifurone standard (g/kg)

**Repeatability  $r$**  = 17.75 g/kg at 960 g/kg active ingredient content

**Reproducibility  $R$**  = 21.74 g/kg at 960 g/kg active ingredient content

---

**FLUPYRADIFURONE ANY OTHER LIQUID  
987/AL/M/-**

**1 Sampling.** Take at least 500 mL\*. Shake the sample well prior to weighing.

\*(for this trial less amount is provided)

**2 Identity tests.**

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

**3 Flupyradifurone.**

Same approach as for Flupyradifurone 987/TC/M/-

**3.1 Determination of Flupyradifurone by reversed phase HPLC**

As for Flupyradifurone 987/TC/M/- except

Disposable PTFE syringe filter compatible with organic solvents and a 0.45 µm pore diameter or centrifuge.

*(c) Calibration solution.* Weigh in duplicate (to the nearest 0.1 mg) approximately 50 mg (s mg) of the Flupyradifurone reference standard into separate volumetric flasks (100 mL). Add 50 mL acetonitrile and place the flasks in an ultrasonic bath for 15 min. Make up the flasks with purified water to just below the calibration mark and allow to cool to ambient temperature. Fill to the mark with purified water and mix thoroughly. Transfer 8 mL of this solution into a separate volumetric flask (100 mL) and make up the flask with acetonitrile / purified water 50 / 50 %v/v to just the calibration mark and allow to cool to ambient temperature. Fill to the mark with acetonitrile / purified water 50 / 50 %v/v and mix thoroughly (Calibration solutions C3, C4) (Fig. 3).

*(d) Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample (*w* mg) to contain about approximately 2 mg (*w* mg) of Flupyradifurone into separate volumetric flasks (50 mL)\*. Add 20 mL acetonitrile and place the flasks in an ultrasonic bath for 15 min. Make up the flasks with acetonitrile to just below the calibration mark and allow to cool to ambient temperature. Fill to the mark with acetonitrile and mix thoroughly. Filter the sample solution through a disposable filter or centrifuge the sample solution (Sample solutions S3, S4) (Fig. 4).

\* This corresponds to a concentration of 4 mg in 100 mL

**Repeatability r** = 0.00180 g/kg at 0.08 g/kg active ingredient content

**Reproducibility R** = 0.00742 g/kg at 0.08 g/kg active ingredient content

---

**FLUPYRADIFURONE EMULSIFIABLE CONCENTRATE**  
**987/EC/M/-**

**1 Sampling.** Take at least 500 mL\*. Shake the sample well prior to weighing.

\*(for this trial less amount is provided)

**2 Identity tests.**

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

**3 Flupyradifurone.**

Same approach as for Flupyradifurone 987/TC/M/-

**3.1 Determination of Flupyradifurone by reversed phase HPLC**

As for Flupyradifurone 987/TC/M/- except

Disposable PTFE syringe filter compatible with organic solvents and a 0.45 µm pore diameter or centrifuge.

*(d) Sample preparation.* Prepare sample solutions in duplicate for each sample. Weigh (to the nearest 0.1 mg) sufficient sample (*w* mg) to contain about approximately 50 mg (*w* mg) of Flupyradifurone into separate volumetric flasks (100 mL). Dissolve in 20 mL purified water and add 50 mL acetonitrile and place the flasks in an ultrasonic bath for 15 min. Make up the flasks with purified water to just below the calibration mark and allow to cool to ambient temperature. Fill to the mark with purified water and mix thoroughly. Filter the sample solution through a disposable filter or centrifuge the sample solution (Sample solutions S5, S6) (Fig. 5).

**Repeatability r** = 0.94 g/kg at 64.7 g/kg active ingredient content

**Reproducibility R** = 1.32 g/kg at 64.7 g/kg active ingredient content

---

**FLUPYRADIFURONE EMULSION, OIL IN WATER  
987/EW/M/-**

**1 Sampling.** Take at least 500 mL\*. Shake the sample well prior to weighing.

\*(for this trial less amount is provided)

**2 Identity tests.**

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

**3 Flupyradifurone.**

Same approach as for Flupyradifurone 987/TC/M/-

**3.1 Determination of Flupyradifurone by reversed phase HPLC**

Same approach as for Flupyradifurone 987/EC/M/-

(Sample solutions S7, S8) (Fig. 6)

**Repeatability r** = 0.94 g/kg at 25 g/kg active ingredient content

**Reproducibility R** = 2.34 g/kg at 25 g/kg active ingredient content

**FLUPYRADIFURONE FLOWABLE CONCENTRATE FOR SEED TREATMENT  
987/FS/M/-**

**1 Sampling.** Take at least 500 mL\*. Shake the sample well prior to weighing.

\*(for this trial less amount is provided)

**2 Identity tests.**

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

**3 Flupyradifurone.**

Same approach as for Flupyradifurone 987/TC/M/-

**3.1 Determination of Flupyradifurone by reversed phase HPLC**

Same approach as for Flupyradifurone 987/EC/M/-

(Sample solutions S9, S10) (Fig. 7)

**Repeatability r** = 8.70 g/kg at 407 g/kg active ingredient content

**Reproducibility R** = 11.71 g/kg at 407 g/kg active ingredient content

## 4 Suspendingibility

REAGENTS AND APPARATUS as for MT 184

### PROCEDURE

(a) *Preparation of suspension and determination of sedimentation. MT 184.*

(b) Determination of Flupyradifurone in the bottom 25 ml of suspension: After removal of the top 225 ml of suspension transfer approx. 250 mg\* out of the remaining and well homogenized 25 ml to a 100 ml volumetric flask and add 20 ml purified water and 50 ml acetonitrile. To dissolve the substance, place the flask in an ultrasonic bath for 15 min and make up the flask with purified water to just below the calibration mark. Wait until the temperature has stabilized, then fill the flask up to the calibration mark with purified water. Determine the mass of flupyradifurone (Q g) by 987/TC/M/-).

\* This is calculated for a 50% suspension, in case of other concentrations the volume has to be adjusted but it has to be ensured to work within the linear range.

## FLUPYRADIFURONE SUSPENSION CONCENTRATE 987/SL/M/-

**1 Sampling.** Take at least 500 mL\*. Shake the sample well prior to weighing.

\*(for this trial less amount is provided)

### 2 Identity tests.

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

### 3 Flupyradifurone.

Same approach as for Flupyradifurone 987/TC/M/-

#### 3.1 Determination of Flupyradifurone by reversed phase HPLC

Same approach as for Flupyradifurone 987/EC/M/-

(Sample solutions S11, S12) (Fig. 8)

**Repeatability r** = 3.20 g/kg at 171 g/kg active ingredient content

**Reproducibility R** = 3.29 g/kg at 171 g/kg active ingredient content



---

**FLUPYRADIFURONE WATER DISPERSABLE GRANULE  
987/WG/M/-**

**1 Sampling.** Take at least 500 g\*. Grind the sample thoroughly in a mortar.

\*(for this trial less amount is provided; please grind the entire sample)

**2 Identity tests.**

**2.1 HPLC.** As for Flupyradifurone 987/TC/M/-

**2.2 UV spectrometry.** As for Flupyradifurone 987/TC/M/-

**3 Flupyradifurone.**

Same approach as for Flupyradifurone 987/TC/M/-

**3.1 Determination of Flupyradifurone by reversed phase HPLC**

Same approach as for Flupyradifurone 987/EC/M/-

(Sample solutions S13, S14) (Fig. 9)

**Repeatability r** = 2.91 g/kg at 120 g/kg active ingredient content

**Reproducibility R** = 3.56 g/kg at 120 g/kg active ingredient content

**4 Suspensibility.** As for Flupyradifurone 987/FS/M/- except

(b) Determination of Flupyradifurone in the bottom 25 ml of suspension. After removal of the top 225 ml of suspension transfer approx. 2700 mg\* out of the remaining and well homogenized 25 ml to a 100 ml volumetric flask and add 20 ml purified water and 50 ml acetonitrile. To dissolve the substance, place the flask in an ultrasonic bath for 15 min and make up the flask with purified water to just below the calibration mark. Wait until the temperature has stabilized, then fill the flask up to the calibration mark with purified water. Determine the mass of flupyradifurone (Q g) by 987/AL/M/-)

\* This is calculated for a 1.2% suspension, in case of other concentrations the volume has to be adjusted but it has to be ensured to work within the linear range.




Fig. 1 UV-Spectrum of Flupyradifurone







Fig. 2 Infrared Spectrum of Flupyradifurone




**Fig. 3 Chromatogram of Flupyradifurone (Standard)**




**Fig. 4 Chromatogram of Flupyradifurone TC**




**Fig. 5 Chromatogram of Flupyradifurone (Standard for AL formulation)**





**Fig. 6 Chromatogram of Flupyradifurone AL**




**Fig. 7 Chromatogram of Flupyradifurone EC**



**Fig. 8 Chromatogram of Flupyradifurone EW**

**Fig. 9 Chromatogram of Flupyradifurone FS****Fig. 10 Chromatogram of Flupyradifurone SL**



**Fig. 11 Chromatogram of Flupyradifurone WG**