

# PYRAOXYSTROBIN

## **Small Scale Collaborative Trial on the Determination of Pyraoxystrobin in Technical Material and Formulations by HPLC**

By

Yu liang


Shenyang Research Institute  
of Chemical Industry  
No.8, Shen Liao Dong Road, Tie Xi district,  
Shenyang City, P.R. of China  
110021

<b>1. INTRODUCTION.....</b>	<b>3</b>
1.1. SCOPE.....	4
1.2. SAMPLES .....	4
1.3. PARTICIPANTS .....	4
<b>2. ANALYTICAL METHODS .....</b>	<b>5</b>
2.1. OUTLINE OF METHODS .....	5
2.2. PROGRAM OF WORK .....	5
<b>3. REMARKS OF PARTICIPANTS.....</b>	<b>6</b>
3.1. ANALYTICAL CONDITIONS .....	6
3.2. REMARKS .....	6
<b>4. RESULTS AND DISCUSSION .....</b>	<b>7</b>
<b>5. CONCLUSION.....</b>	<b>7</b>

## 1. INTRODUCTION

### Pyraoxystrobin

964



*ISO Common name* Pyraoxystrobin

*CAS No* [862588-11-2]

*Chemical name:*

(*E*)-methyl

2-(2-((3-(4-chlorophenyl)-1-methyl-1*H*-pyrazol-5-yloxy)methyl)phenyl)-  
3-methoxyacrylate

*Empirical formula* C<sub>22</sub>H<sub>21</sub>ClN<sub>2</sub>O<sub>4</sub>

*RMM* 412.87

*m.p.* 129.6°C

*v.p.* 1.22E-008Pa

*Solubility* In water(g/L,20°C): 0.03mg/L, high solubility  
in DMF, acetone, ethyl acetate, methanol;  
Low solubility in petro-ether and water.

<i>Description</i>	The pure material is a white, odourless solid
<i>Stability</i>	Not stable in alcohol to heat
<i>Formulations</i>	Suspension concentrates

### **1.1. Scope**

The results of the small scale collaborative study for pyraoxystrobin technical material and formulations are reported.

### **1.2. Samples**

- 1) Pyraoxystrobin technical (TC-1)
- 2) Pyraoxystrobin suspension concentrates (SC-1)
- 3) Pyraoxystrobin suspension concentrates (SC-2)

### **1.3. Participants**

- 1) Ms. Sun menglin Bao ning

SHENYANG SCIENCREAT CHEMICALS  
CO.,LTD.

No.17, Xihe North street, economic and  
technology development zone , ShenYang City,  
P.R.china.

- 2) Mr. Xian shuijin

Nutrichem Laboratory Co.,Ltd

D-1 , Dongsheng Science Park  
66 Xi xiaokou Road, Haidan District, Beijing  
100192, P.R.china

- 3) Mr. Zou weibin

JIANGSU ROTAM CHEMISTRY CO., LTD

No.88 Rotam Rd. ETDZ, Kunshan,

jiangsu,215301 P.R.china.

4) Ms Zhang xuebing                      Pesticides Test Laboratory of Shenyang Research  
Institute of Chemical Industry  
No.8, Shen Liao Dong Road, Tie Xi district,  
ShenYang City, P.R.china

5) Ms. Li xiujie                              Shenyang Research Institute of Chemical Industry  
No.8, Shen Liao Dong Road, Tie Xi district,  
Shenyang City, P.R. of China  
110021

## **2. ANALYTICAL METHODS**

### **2.1. Outline of Methods**

Pyraoxystrobin in the test samples is determined by HPLC using uv detection.

Pyraoxystrobin Analytical standard, Batch No.: 20110927, Purity 99.7%;

### **2.2. program of work**

We requested the collaborators to:

- 1) conduct duplicate determinations on two different days for each of the three samples;
- 2) inject each sample solution in duplicate and calculate the mean value;
- 3) describe operating conditions in detail;
- 4) recommended procedure if the method is deviate .

### 3. Remarks of participants

#### 3.1. Analytical conditions

Lab	Liquid chromatograph integrator	Column	Mobile phase	Flow rate (ml/min)	Column temp(°C)
1	Waters Empower2	Eclipse XDB-C18 150mm× 4.6mm	Acetonitrile/H <sub>2</sub> O 60:40	1.3	ambient
2	Agilent HPLC ChemStation	Zorbax XDB-C18 150mm x 4.6 mm	Acetonitrile/H <sub>2</sub> O 60:40	1.3	20
3	Agilent HPLC ChemStation	Agilent XDB-C18 150mm× 4.6mm	Acetonitrile/H <sub>2</sub> O 60:40	1.3	ambient
4	Agilent HPLC ChemStation	Agilent Zorbax C18 150mm x 4.6 mm	Acetonitrile/H <sub>2</sub> O 60:40	1.5	ambient
5	Waters Empower2	SB-C18 150mm× 4.6mm	Acetonitrile/H <sub>2</sub> O 60:40	1.3	ambient

#### 3.2. Remarks

- Lab 1

The column temperature had changed more than 5°C.

- Lab 4

The sample SC-2 was not dissolved completely when it Ultrasonic less than 15min.

#### 4. Results and Discussion

Five data sheets were obtained from the five participants. Summary and detailed statistical evaluations are shown in Tables 1, 2-1 and 2-2. The statistical evaluations were carried out according to ISO 5725.

There is no stragglers and outliers in the small scale collaborative trial on Pyraoxystrobin .

#### 5. CONCLUSION

For all samples, the values of  $RSD_R$ (reproducibility relative standard deviation) were smaller than those calculated by Horwitz' s equation. The proposed method is considered appropriate for the determination of pyraoxystrobin in technical and suspension concentrates product.

We proposes proceeding to a large scale collaborative study.

Table 1 Summary of Statistical Evaluation of Pyraoxystrobin Small Scale Collaborative Study

	<i>TC-1</i>	<i>SC-1</i>	<i>SC-2</i>
<i>Average(g/kg)</i>	973.78	203.56	199.06
<i>Number of labs.</i>	5	5	5
<i>Repeatability standard deviation(<math>S_r</math>)</i>	2.145	1.200	0.951

<b><i>“Pure ” between laboratory standard variation(<math>S_L</math>)</i></b>	4.785	0.857	2.210
<b><i>Reproducibility standard deviation(<math>S_R</math>)</i></b>	5.244	1.475	2.405
<b><i>Repeatability(<math>r</math>)</i></b>	6.006	3.360	2.663
<b><i>Reproducibility(<math>R</math>)</i></b>	14.683	4.130	6.734
<b><i><math>RSD_r</math></i></b>	0.220	0.590	0.478
<b><i><math>RSD_R</math></i></b>	0.539	0.725	1.208
<b><i>Horwitz’s value</i></b>	2.008	2.542	2.551

Table 2-1 Statistical Evaluation (Technical -1)

<i>Lab.</i>	<i>Analytical data</i>					
	<i>(n=4)</i>		<i>Y<sub>i</sub></i>	<i>(Y<sub>i</sub>)<sup>2</sup></i>	<i>S<sub>i</sub></i>	<i>Si<sup>2</sup></i>
1 Day1	979.7	980.4	980.32	961027.30	0.592	0.350
	980.1	981.1				
2 Day1	975.0	980.3	977.50	955506.25	3.492	12.194
	974.0	980.7				
3 Day1	972.0	972.7	972.08	944939.53	0.436	0.190
	971.7	971.9				
4 Day1	968.6	973.4	969.80	940512.04	2.903	8.427
	966.6	970.6				
5 Day1	970.8	967.7	969.22	939387.41	1.358	1.844



Day2	969.8	968.6
<b>S1 SUM Yi=</b>	4868.92	
<b>S2 SUM Yi<sup>2</sup>=</b>	4741372.50	
<b>S3 SUM Si<sup>2</sup></b>	23.005	

1) Cochran's test (p=5,n=4)

$$C = S_i^2 \max / S3 = 0.530 < 0.598 \text{ (p=5,n=4, 5\%)}$$

2) Grubb's test (p=5,n=4)

$$Y_i(\min) = 969.22 \quad Y_i(\max) = 980.32 \quad Y = S1/p \quad 973.78$$

$$S = 4.902$$

$$Y - Y_i(\min) = 4.56$$

$$Y_i(\max) - Y = 6.54$$

$$\text{Lower} = [Y - Y_i(\min)] / S = 0.930 < 1.67 \text{ (p=5,n=4, 5\%)}$$

$$\text{Upper} = [Y_i(\max) - Y] / S = 1.334 < 1.67 \text{ (p=5,n=4, 5\%)}$$

3) Calculation of r and R

$$\text{Mean; } Y = S1/p = 973.78$$

$$S_r^2 = S3 / p = 4.601 \quad S_r = 2.145$$

$$S_L^2 = [(pS2 - S1^2) / p(p-1)] - (S_r^2/n) = 22.900 \quad S_L = 4.785$$

$$S_R^2 = S_r^2 + S_L^2 = 27.501 \quad S_R = 5.244$$

$$r = 2.8 \times S_r = 6.006$$

$$R = 2.8 \times S_R = 14.683$$

$$RSDr = (S_r / \text{mean}) \times 100 = 0.220$$

$$RSDR = (S_R / \text{mean}) \times 100 = 0.539$$

$$\text{Horwitz's Value} = 2^{[1 - 0.5 \times \log(Y/1000)]} = 2.008$$

$$\text{RSDR} < 2.008 \quad (\text{Horwitz's Value})$$

Table 2-2 Statistical Evaluation (SC -1)

<i>Lab.</i>	<i>Analytical data (n=4)</i>		<i>Yi</i>	<i>(Yi)<sup>2</sup></i>	<i>Si</i>	<i>Si<sup>2</sup></i>
1 Day1	205.7	204.1	203.40	41371.56	1.937	3.752
	202.6	201.2				
2 Day1	204.7	204.9	205.20	42107.04	1.543	2.381
	207.4	203.8				
3 Day1	203.8	204.2	203.60	41452.96	0.479	0.229
	203.4	203.1				
4 Day1	203.3	202.9	203.30	41330.89	0.432	0.187
	203.9	203.1				
5 Day1	202.7	203.2	202.30	40925.29	0.810	0.656
	201.5	201.7				
<b>S1 SUM Yi=</b>			1017.80			
<b>S2 SUM Yi<sup>2</sup>=</b>			207187.74			
<b>S3 SUM Si<sup>2</sup></b>						7.205

1) Cochran's test (p=5, n=4)

$$C = Si^2_{\text{max}} / S3 = 0.521 < 0.598 \quad (p=5, n=4, 5\%)$$

2) Grubb's test (p=5, n=4)

$$Y_i(\min) = 202.30 \quad Y_i(\max) = 205.20 \quad Y = S1/p \quad 203.56$$

$$S=1.045$$

$$Y - Y_i(\min) = 1.26$$

$$Y_i(\max) - Y = 1.64$$

$$\text{Lower} = [Y - Y_i(\min)] / S = 1.206 < 1.67 \text{ (p=5, n=4, 5\%)}$$

$$\text{Upper} = [Y_i(\max) - Y] / S = 1.569 < 1.67 \text{ (p=5, n=4, 5\%)}$$

### 3) Calculation of r and R

$$\text{Mean; } Y = S1/p = 203.56$$

$$S_r^2 = S3 / p = 1.441 \quad S_r = 1.200$$

$$S_L^2 = [(pS2 - S1^2)/p(p-1)] - (S_r^2/n) = 0.735 \quad S_L = 0.857$$

$$S_R^2 = S_r^2 + S_L^2 = 2.176 \quad S_R = 1.475$$

$$r = 2.8 \times S_r = 3.360$$

$$R = 2.8 \times S_R = 4.130$$

$$RSD_r = (S_r / \text{mean}) \times 100 = 0.590$$

$$RSD_R = (S_R / \text{mean}) \times 100 = 0.725$$

$$\text{Horwitz's Value} = 2 [1 - 0.5 \times \log(Y / 1000)] = 2.542$$

$$RSD_R < 2.542 \quad (\text{Horwitz's Value})$$

Table 2-3 Statistical Evaluation (SC -2)

<i>Lab.</i>	<i>Analytical data (n=4)</i>		<i>Y<sub>i</sub></i>	<i>(Y<sub>i</sub>)<sup>2</sup></i>	<i>S<sub>i</sub></i>	<i>S<sub>i</sub><sup>2</sup></i>
1 Day1	197.4	196.9	196.02	38423.84	1.325	1.756
Day2	194.7	195.1				
2 Day1	196.0	199.5	197.55	39026.00	1.515	2.295

Day2	196.8	197.9				
3 Day1	202.1	201.5				
			201.62	40650.62	0.320	0.102
Day2	201.5	201.4				
4 Day1	199.5	198.9				
			199.60	39840.16	0.577	0.333
Day2	199.7	200.3				
5 Day1	200.7	200.6				
			200.50	40200.25	0.183	0.033
Day2	200.3	200.4				
<b>S1 SUM Yi=</b>			995.29			
<b>S2 SUM Yi<sup>2</sup>=</b>				198140.87		
<b>S3 SUM Si<sup>2</sup></b>						4.519

1) Cochran's test (p=5, n=4)

$$C = S_i^2 \max / S_3 = 0.508 < 0.598 \text{ (p=5, n=4, 5\%)}$$

2) Grubb's test (p=5, n=4)

$$Y_i(\min) = 196.02 \quad Y_i(\max) = 201.62 \quad Y = S_1/p = 199.06$$

$$S = 2.261$$

$$Y - Y_i(\min) = 3.04$$

$$Y_i(\max) - Y = 2.56$$

$$\text{Lower} = [Y - Y_i(\min)] / S = 1.345 < 1.67 \text{ (p=5, n=4, 5\%)}$$

$$\text{Upper} = [Y_i(\max) - Y] / S = 1.132 < 1.67 \text{ (p=5, n=4, 5\%)}$$

3) Calculation of r and R

$$\text{Mean; } Y = S_1 / p = 199.06$$

$$S_r^2 = S_3 / p = 0.904$$

$$S_r = 0.951$$

$$SL^2 = [(pS^2 - S_1^2)/p(p-1)] - (Sr^2/n) = 4.882 \quad SL = 2.210$$

$$SR^2 = Sr^2 + SL^2 = 5.786 \quad SR = 2.405$$

$$r = 2.8 \times Sr = 2.663$$


$$R = 2.8 \times SR = 6.734$$

$$RSDr = (Sr / \text{mean}) \times 100 = 0.478$$


$$RSDR = (SR / \text{mean}) \times 100 = 1.208$$

$$\text{Horwitz's Value} = 2 [1 - 0.5 \times \log (Y / 1000)] = 2.551$$


$$RSDR < 2.551 \quad (\text{Horwitz's Value})$$



**Fig.1 Praoxystrobin Technical-1**



**Fig.2 Praoxystrobin SC-1**



**Fig.3 Praoxystrobin SC-2**