

## Handbook 1C

<b>Compound</b>	<b>Method</b>	<b>Handbook</b>	<b>page</b>	<b>Method based on</b>
Azinphos-methyl		1C		1970 Method extension for WP and EC based on the AOAC quantitative IR method in Handbook 1
BHC	[4], method			1977 WHO CIPAC Method based on packed GC on OV-210
Bromacil	[139], method	1C		1986 titration
Bromoxynil octanoate	GC method	1C		1998 packed GC OV 101
Bromoxynil octanoate	GC method	1C		2004 AOAC packed GC SP-2100
Bupirimate	[261], method	1C		2007 GC packed column SE30
Chinomethionat	[172], method	1C		2019 HPLC NP
Chlorthal-dimethyl	[328], method	1C		2032 AOAC packed GC on UC98, includes determination of hexachlorobenzene as impurity
Chlorthal-dimethyl		1C		2034 quantitative IR method
Cypermethrin	[332], method, GC	1C		2047 packed GC, OV101
2,4-D-esters	-esters titration method	1C		2068 titration
1,2-Dibromoethane	[125], method	1C		2115 AOAC, packed GC, tricresyl phosphate
Dichlorprop	GC	1C		2087 packed GC, Apiezon L
Dichlorprop-salts	HPLC	1C		2093 RP HPLC
Diclofop-methyl	[358], method	1C		2096 packed GC after silylation on OV101
Disulfoton	[152], method	1C		2101 AOAC, packed GC, SE30
Dithianon	[153], method	1C		2105 colorimetric after reaction with borate
Endosulfan	[89], method	1C		2111 packed GC, OV210
Fenitrothion	[35], method	1C		2117 WHO, packed GC, PPE-6R
Fenoprop	GC	1C		2122 packed GC, Apiezon L
Fensulfotion	[303], method	1C		2126 method, HPLC
Formaldehyde	[156], method	1C		2130 titration after oxidation to formic acid
MCPA	[2], titration method	1C		2137
MCPA	GC	1C		2138 packed GC, Apiezon L
MCPA-salts	GC	1C		2141 packed GC, Apiezon L
MCPA-esters	- esters, titration method	1C		2145

## Handbook 1C

<b>Compound</b>	<b>Method</b>	<b>Handbook</b>	<b>page</b>	<b>Method based on</b>
MCPA-esters	GC	1C	2145	packed GC, Apiezon L
MCPA + 2,3,6-TBA	[2 + 49]	1C	2148	packed GC, Apiezon L
MCPB	[50], titration method	1C	2153	titration
MCPB	GC	1C	2154	packed GC, Apiezon L
MCPB-salts	- salts, titration method	1C	2154	titration
MCPB-salts	GC	1C	2155	packed GC, Apiezon L
Mecoprop	[51] GC	1C	2158	packed GC, Apiezon L
Mecoprop-salts	GC	1C	2163	packed GC, Apiezon L
MGK 264 ,	[345]	1C	2166	packed GC, OV-1
Parathion	[10]	1C	2169	extension to CS formulations, AOAC, packed GC, OV17
Parathion-methyl	[487]	1C	2169	extension to CS formulations, AOAC, packed GC, OV17
Piperonyl butoxide	[33], method	1C	2190	AOAC, packed GC, OV101
Potassium cyanate [197]	[197], method	1C	2200	AOAC, gravimetric method
Prometon	[175], method	1C	2202	titration with perchloric acid in glacial acetic acid
Prometon		1C	2203	AOAC, packed GC, Carbowax
Pyrazophos	[350], method	1C	2206	NP HPLC
Pyrethrum	[32], method	1C	2209	AOAC, packed GC, OV101
Pyrethrum + piperonylbutoxide	[32 + 33]	1C	2209	AOAC, packed GC, OV101
Quintozene	[78], method	1C	2213	AOAC, packed GC, SE30
2,4,5-T	[6], method	1C	2220	titration
2,4,5-T	GC	1C	2221	packed GC, Apiezon L
2,3,6-TBA	[49], method	1C	2225	titration
2,3,6-TBA		1C	2227	packed GC
2,3,6-TBA Na-salt	Na-salt titration method	1C	2228	titration

Handbook 1C

<b>Compound</b>	<b>Method</b>	<b>Handbook</b>	<b>page</b>	<b>Method based on</b>
2,3,6-TBA Na-salt	GLC	1C	2229	packed GC
Thallosulfate	181], method	1C	2234	AOAC, gravimetry